



BARRIER®
Staff Clothing
Product Guide
2022

BARRIER® Staff clothing

The protection you need. The comfort you want.

Top 4 reasons



Extra comfort



Infection prevention⁶



Cost effective



Environmentally conscious

Unparalleled comfort with our exclusive Unisoft® material

- Better air permeability and wicks moisture away from skin faster than reusable scrubs¹
- Less linting particles
- Lightweight
- Soft feel of cotton
- Anti-static treated
- Quiet
- Breathable
- Not made with natural rubber latex

Risk reduction for you and your patients

80% of bacteria come from the air.² You shed thousands of skin squames every minute,³ each potentially carrying bacteria.


5% of patients risk a surgical site infection⁴

SSI's are causing trauma and unnecessary suffering for patients, and costing up to **\$10 billion** annually⁵



3X
— and —
4X

BARRIER's exclusive Unisoft material is spunbonded rather than woven, which results in a reduced pore size, providing up to three times greater bacterial protection and up to four times greater fluid protection than leading reusable scrubs.⁶

Reusable Woven Fabric



BARRIER Unisoft material



BARRIER meets AORN guidelines for surgical attire



Reduce overall scrub costs

Nearly **90%** of all linen, like scrubs, used in U.S. hospitals does not reach its useful life.⁷



BARRIER® Extra Comfort single-use scrubs and jackets can help reduce your overall scrub costs by eliminating processing costs, replacement costs and laundering costs.



Conscious of the environment

The Swedish Environmental Research Institute performed an assessment of the environmental impacts of BARRIER Extra Comfort scrubs and cotton/polyester reusable scrubs. Material production, scrub suit production, distribution, laundry, waste treatment and emissions were all considered.

BARRIER scrubs have a smaller environmental impact compared to reusable scrubs⁸

- Lower CO₂ emissions (climate change)
- Less inorganic matter introduced into our bodies of water (eutrophication)
- Less pollutants into the air (ground level ozone)

BARRIER Scrub Suits can be better for the environment.




Table of Contents

- Scrub Suits 3
- Long-sleeved Jackets 5
- Protective Eyewear 6
- Surgical Headwear 7
- Standard Surgical Face Masks 8
- Special Surgical Face Masks 9
- Extra Protection Surgical Face Masks 10

The information provided herein is not to be construed as the practice of medicine or substituted for the independent medical judgment of a patient's treating clinician. This information, including but not limited to suggestions for product selection and suggested use is based on generalizations and does not consider the unique characteristics of an individual patient. Each patient's clinician shall remain solely responsible for assessing patient, determining the appropriate treatment and managing treatment. For additional information, please refer to the applicable product insert or contact Mölnlycke Health Care at 1-800-843-8497.

BARRIER® Extra Comfort Scrub Suits

Single-use scrubs offer a new, fresh start every day



Help drive your infection control effort. BARRIER scrubs are produced with a low-linting material in a controlled clean environment so you will always know that they are hygienically clean. With traditional reusable scrubs you can never really know for sure.

Fashionable design combined with a professional look.

Durability to last all day
Strong 100% polypropylene material with reinforced stitching to minimize tearing, even in high stress areas.

Always fresh and clean.
No wear and tear.
No stains or odors.
No need to worry.

Available in two colors*

- Green
- Blue

Large size range*

from XS through 4XL

Separately packed

shirts and pants allows for individual selection of sizes



* Not all colors and styles are available in every size, see ordering information below

Ordering Information

Extra Comfort Scrub Shirts

Description	XS	S	M	L	XL	2XL	3XL	4XL
Scrub Shirt, Green		18610	18620	18630	18640	18650	18660	18670
Scrub Shirt, Blue	21600	21610	21620	21630	21640	21650	21660	21670

48 pieces per case. 4 inner packs, each containing 12 pieces

Extra Comfort Scrub Pants

Description	XS	S	M	L	XL	2XL	3XL	4XL
Drawstring Pants, Green		18710	18720	18730	18740	18750	18760	18770
Drawstring Pants, Blue	21700	21710	21720	21730	21740	21750	21760	21770
Elastic Waist Pants, Green		18920	18930	18940	18950	18960		

48 pieces per case. 4 inner packs, each containing 12 pieces

BARRIER® Extra Comfort Long-Sleeved Jacket

Single-use jackets offer a new, fresh start every day



The AORN Guidelines for Surgical Attire recommends that all non-scrubbed personnel wear a freshly laundered or single-use long-sleeved jacket snapped closed with cuffs down to the wrists.³

Long-sleeved attire helps decrease the risk of airborne transmission of bacteria from the shedding of skin squames.³

One-day use only jackets reduces the chance for the daily accumulation of harmful bacteria.

Ordering Information

Description	S	M	L	XL	2XL	3XL
Long-sleeved Jacket, Green	18001	18010	18020	18030	18040	18050
Long-sleeved Jacket, Blue	28000	28010	28020	28030	28040	28050
XP Long-sleeved Jacket, Green	18200	18210	18220	18230	18240	18250

48 pieces per case. 4 inner packs, each containing 12 pieces

BARRIER® Protective Eyewear

Protective glasses



Features

- Perfect for use in the OR
- Lightweight and comfortable
- Designed for maximum peripheral vision without distortion
- Solid side shields per OSHA guidelines
- Anti-fog lens



Ordering Information

Description	Product Code	Box	Case
Protective Glasses (one size)	1702	10	30


BARRIER® Headwear

Surgical caps and hoods

A B Sheer bouffant nurse's cap



C Operating cap



D Surgical cap w/ tie strings



E Operating cap



Features

- Breathable for better comfort
- Unique elastic qualities combine flexibility and softness
- Low-linting
- Strong and durable material helps contain shedding from hair particles
- Not made with natural rubber latex

Ordering Information

Description	M	L	Box	Case
A Sheer Bouffant Nurse's Cap, Green	621725	621825	150	600
B Sheer Bouffant Nurse's Cap, Blue	621715	621815	150	600
Description	Product code	Box	Case	
C Operating Cap (one size)	620500	100	500	
D Surgical Cap w/ Tie Strings (one size)	621301	100	500	
E Surgical Hood (one size)	620105	100	500	

BARRIER® Surgical Face Masks

Standard surgical face masks

A Surgical mask w/ ear loop



B Surgical Mask



C Surgical Mask



Features

- Malleable nose wire
- Pleated, 3-layer design
- Not made with natural rubber latex
- Tie-band or ear loop
- Extra long tie bands

Specifications

- Standard EN 14683 Type II
- Class I Flammability rating
- Bacterial Filtration Efficiency (BFE) > 98%
- Low breathing resistance

Ordering Information

Description	Product code	Box	Case
A Surgical Mask w/ Ear Loop	42281	50	500
B Surgical Mask, Blue	42301	60	600
C Surgical Mask, Green	42391	60	600

BARRIER® Surgical Face Masks

Special surgical face masks

A Anti-fog surgical mask



Features

- Sure-fit nose clip
- Anti-fog
- Not made with natural rubber latex
- Pleated 3-layer design or duckbill
- Soft cellulose inner facing layer
- Tie band or ear-loop
- Extra long tie bands

B Anti-fog surgical mask



Specifications

- Standard EN 14683 Type II
- Class I Flammability rating
- Bacterial Filtration Efficiency (BFE) > 98%
- Low breathing resistance


Ordering Information

	Description	Product code	Box	Case
A	Anti-fog Surgical Mask, Blue	42311	60	600
B	Anti-fog Surgical Mask, Green	42381	60	600

BARRIER® Surgical Face Masks

Extra protection surgical face masks

A Surgical mask w/ ear loop



C D Surgical mask w/ visor



E Surgical Mask



F Surgical Mask



Features


- Sure-fit nose clip
- Anti-fog
- Not made with natural rubber latex
- Soft cellulose inner facing layer
- Tie band or ear-loop
- Extra long tie bands
- 4-layer, splash resistant design

Specifications

- Standard EN 14683 Type IIR
- Class I Flammability rating
- Bacterial Filtration Efficiency (BFE) > 98%
- Fluid resistant 160 mm Hg

Ordering Information

	Description	Product code	Box	Case
A	Anti-fog, Splash-resistant Surgical Mask w/ Ear Loop, Blue	42241	50	500
B	Anti-fog, Splash-resistant Surgical Mask w/ Visor	42321	50	200
C	Anti-fog, Splash-resistant Surgical Mask w/ Ear Loop & Visor	42671	50	200
D	Anti-fog, Splash-resistant Surgical Mask, Blue	42341	50	500
E	Anti-fog, Splash-resistant Surgical Mask, Lavender	42361	50	500



For further information
on BARRIER® products
or to request samples,
please visit
www.molnlycke.us
or email us at
info.us@molnlycke.com

References: **1.** Molnlycke Laboratory Test, PD-528943.cc **2.** Howarth FW. Prevention of airborne infection during surgery. *Lancet*. 1985; 1(8425): 386-8. **3.** Association of periOperative Registered Nurses (AORN) Guidelines for Perioperative Practice, Guideline for Surgical Attire (2017). **4.** Cheadle WG, Risk factors for surgical site infection. *Surg Infect (Larchmt)*. 2006;7 Suppl 1: S7-11. **5.** Loyola University Health System. Surgical site infections are the most common and costly of hospital infections: Guidelines for preventing surgical site infections are updated. *ScienceDaily*, 19 January 2017. www.sciencedaily.com/releases/2017/01/170119161551.htm **6.** Molnlycke Health Care. Data on File, 2014. **7.** Stern, David. Linen loss: the \$850 million elephant in your rooms. *McKnight's Long-Term Care News*. Sept. 10 2014. <http://www.mcknights.com/marketplace/linen-loss-the-840-million-elephant-in-your-rooms/article/367641/> **8.** Eriksson E & Westerdahl J. Life cycle assessment of single-use and reusable hospital staff clothing on the Spanish market. Gothenburg, 2 September, 2011. IVL Swedish Environmental Research Institute.

BARRIER®

The Molnlycke, BARRIER trademarks, names and logo types are registered globally to one or more of the Molnlycke Health Care Group of Companies. Distributed by Molnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2022 Molnlycke Health Care AB. All rights reserved. 1-800-843-8497 MHC-2021-79168

Lit# 1-0360 09.20


Molnlycke®