

PRODUCT DATA SHEET

MEFIX® PHARMACY ASSORTMENT

Product name	Mefix						
Product description	Self adhesive fabric						
Intended use	Self-adhesive fabric for fixation of dressings, tubes and cannulas etc.						
Product Codes	Art no*	Description	Size (cm x m)	Size (inch)	Pack count con/ret/trp	Volume Outer dimensions (mm)	Volume Outer dimensions (mm)
	310580-05	Mefix 5cmx2.5m ph Fr	5 x 2.5	2x100	1/1/30	56x59x60	295x195x157
	310586-05	Mefix 5cmx5m ph FR	5 x 5	2x200	1/1/30	56x59x60	295x195x157
	310540-00*	Mefix 5cmx5m ph ES	5 x 5	2x200	1/1/24	62x62x62	263x195x157
	310584-05	Mefix 5cmx10m ph FR	5 x 10	2x400	1/1/40	56x74x70	395x236x157
	311080-05	Mefix 10cmx2.5m Ph FR	10 x 2.5	4x100	1/1/30	106x59x57	295x195x246
	311086-05	Mefix 10cmx5m Ph FR	10 x 5	4x200	1/1/30	106x59x57	295x195x246
	311084-05	Mefix 10cmx10m ph FR	10 x 10	4x400	1/1/20	106x74x70	395x236x157
	311040-00*	Mefix 10cmx10m ph ES	10 x 10	4x400	1/1/30	112x75x75	395x236x246
	311580-05	Mefix 15cmx2.5m ph FR	15x2.5	6x100	1/1/20	156x59x57	295x195x246
	311586-05	Mefix 15cmx5m ph FR	15x5	6x200	1/1/20	156x59x57	295x195x246
	*) For the Spanish assortment the Mefix product comes on a plastic core						
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden						
Country of Origin	Finland						

Product classification according to MDD	Class I			
CE-marked	Yes			
Sterilization method	N.A. Mefix is a non sterile product			
Shelf life	5 years			
Products material	Non woven	Spun laced non woven material that consist of polyester fibres		
	Adhesive	Water based poly acrylic adhesive		
	Release	Paper release that consist of siliconised paper		
	Core Only available on 310540-00, 311040-00	Plastic; 100% polystyrene		
Environmental classification	Mefix product	Green		
	According to the ISO 14001 all raw materials are classified from an environmental point of view.			
Other info of special interest	LATEX: Mefix Products and packaging are latex free.			
Physical characteristics	Property	Value	Reference method	
	Water vapour permeability	>6500 g/m ² /24h	ASTM D 6701-01	
Packaging	The purpose of the dispenser box and transport box is to contain and protect material from dust, moisture and damage during transport, handling and storage. Unit package is a sterile barrier for the product; sterility is guaranteed unless inner package is damaged or opened prior to use.			
Storage Condition	Mefix articles should be stored under normal conditions (8-40°C, 46-104°F, 30-70% humidity)			
GS1 Codes	Art No	CON	RET	TRP
	310580-05	-	7332551829643	7332551829636
	310586-05	-	7332551829759	7332551829742
	310540-00	-	7332430754929	7332430754912
	310584-05	-	7332551829797	7332551829780
	311080-05	-	7332551829858	7332551829841
	311084-05	-	7332551829940	7332551829933
	311040-00	-	7332430754967	7332430754950
	311086-05	-	7332551829919	7332551829902
	311580-05	-	7332551830038	7332551830021
	311586-05	-	7332551830090	7332551830083