



# Tubifast® 2-Way Stretch®: What size to choose?


## RED LINE




## GREEN LINE




## BLUE LINE



## YELLOW LINE



## PURPLE LINE



Limb	Limb circumference	Product line	Product length
Small limbs	9 – 18 cm   3.5 – 7 in	<b>RED LINE</b>	1 m – 10 m
Small & medium limbs	14 – 24 cm   5.5 – 9.5 in	<b>GREEN LINE</b>	1 m – 3 m – 5 m – 10 m
Large limbs	24 – 40 cm   9.5 – 16 in	<b>BLUE LINE</b>	1 m – 3 m – 5 m – 10 m
Extra-large limbs, heads, children's trunks	35 – 64 cm   14 – 25 in	<b>YELLOW LINE</b>	1 m – 3 m – 5 m – 10 m
Adult trunks	64 – 130 cm   25 – 51 in	<b>PURPLE LINE</b>	1 m – 5 m – 10 m

Find Tubifast® 2-Way Stretch® application videos at <https://www.molnlycke.ca/education/wound-areas/wound-healing/cutting-guide-for-dressings/>

Find out more at [www.molnlycke.ca](http://www.molnlycke.ca)

Mölnlycke Health Care Canada. 2010 Winston Park Drive, Suite 100 Oakville, ON L6H 5R7. Phone: 905-829-1502  
The Mölnlycke Health Care, Tubifast and TwoWay Stretch trademarks, names and respective logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies © 2023 Mölnlycke Health Care AB. All rights reserved.

