

PRODUCT DATA SHEET
MEPITAC®

Product name	Mepitac				
Product description	Mepitac is a non-sterile soft silicone tape				
Intended use	Mepitac is designed for fixation of medical devices such as drains, tubes, probes, electrodes, IV cannulae and dressings. Mepitac offers gentle skin protection when used under devices e.g. tubes.				
Product Codes	Art.No	Product Name	Size	Pack count CON/RET/TRP	
	298300-00	Mepitac 2cm x 3m	2cm x 3m	1/1/12	
	298370-00	Mepitac 2cm x 3m OTC Norden	2cm x 3m	0/1/12	
	298400-00	Mepitac 4cm x 1.5m	4cm x 1.5m	1/1/12	
	298488-00	Mepitac 0,8in x 3,3yd / 2cm x 3m, CVS	0,8in x 3,3yd / 2cm x 3m	1/1/12	
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden				
Country of origin	Finland				
Product classification according to MDD	Class I, according to Rule 1				
CE-marked	Yes, Notified Body BSI 0086				
Products material	Backing	Polyurethane film / Polyacrylate adhesive / Polyester/polyamide tricot			
	Wound contact layer	Soft Silicone			
	Release film	Polyethylene film			
Shelf-life	3 years				
Physical /Chemical characteristics	Mepitac is released according to the following standards				
	Property	Result	Reference Method		
	Adhesion	≥ 0,7 N/20mm ≤ 2,4 N/20mm	T-261		
Packaging	The purpose of the shelf container and transport box is to protect material from dust, moisture and damage during transport, handling and storage.				
Storage Condition	N/A				

GS1 Codes	Art No	CON	RET	TRP
	298300-00	7332430470836	7332430470836	7332430409836
	298370-00	-	7323190124855	7323190124848
	298400-00	7332430470843	7332430470843	7332430409867
	298488-00	7323190176557	7323190155170	7323190155163